


# Non-Compliant PTCA Balloon Catheter Compliance Chart

Balloon Length 8 mm – 20 mm



Pressure (bar)	Balloon Diameter (mm)										
	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.50	5.00
8	1.96	2.20	2.45	2.69	2.94	3.18	3.42	3.67	3.90	4.35	4.92
10	1.98	2.23	2.47	2.72	2.97	3.22	3.46	3.71	3.95	4.43	4.96
12	2.00	2.25	2.50	2.75	3.00	3.25	3.50	3.75	4.00	4.50	5.00
14	2.02	2.27	2.53	2.78	3.03	3.29	3.54	3.79	4.05	4.55	5.07
16	2.04	2.30	2.56	2.81	3.07	3.32	3.58	3.83	4.10	4.62	5.14
18	2.07	2.33	2.58	2.84	3.10	3.36	3.63	3.88	4.16	4.69	5.22
20	2.09	2.35	2.61	2.87	3.13	3.40	3.67	3.92	4.21	4.75	5.28
22	2.11	2.38	2.64	2.90	3.17	3.43	3.72	3.96	4.27	4.82	5.35
24	2.14	2.40	2.67	2.94	3.20	3.47	3.76	4.00	4.32	4.89	5.40
26	2.16	2.43	2.70	2.97	3.24	3.51	3.80	4.05	4.38	4.97	5.49
28	2.18	2.46	2.73	3.00	3.27	3.55	3.85	4.09	4.44	5.00	5.57

 **Nominal Balloon Pressure**

 **Rated Burst Pressure**  
 20 bar - Ø 4.50 - 5.00 mm  
 22 bar - Ø 2.00 - 4.00 mm

The nominal in vitro device specifications do not take into account any lesion resistance. Do not exceed RBP.