

# Semi-Compliant PTCA Balloon Catheter Compliance Chart

Balloon Length 10 mm - 30 mm



Pressure (bar)	Balloon Diameter (mm)									
	1.00	1.25	1.50	1.75	2.00	2.25	2.50	3.00	3.50	4.00
2	0.90	1.15	1.40	1.65	1.87	2.15	2.36	2.82	3.30	3.78
4	0.95	1.20	1.45	1.70	1.94	2.20	2.45	2.91	3.40	3.90
6	1.00	1.25	1.50	1.75	2.00	2.25	2.50	3.00	3.50	4.00
8	1.05	1.30	1.55	1.80	2.05	2.35	2.57	3.10	3.63	4.13
10	1.10	1.35	1.60	1.85	2.11	2.40	2.64	3.18	3.74	4.25
12	1.15	1.40	1.65	1.90	2.19	2.50	2.71	3.27	3.85	4.35
14	1.20	1.45	1.70	1.95	2.25	2.55	2.80	3.35	3.95	4.45
16	1.23	1.50	1.75	2.00	2.32	2.60	2.85	3.45	4.05	4.58
18	1.25	1.55	1.80	2.05	2.38	2.68	2.93	3.55	4.15	4.68

 Nominal Balloon Pressure

 Rated Burst Pressure

The nominal in vitro device specifications do not take into account any lesion resistance. Do not exceed RBP.